

ECA Check® Plus Easygel®

Specification Data

Selective For:	<i>E. coli</i> , total coliforms, <i>Salmonella spp.</i> and <i>Aeromonas spp.</i>
Description:	ECA Check® Plus is a clear medium containing a chromogenic substrate mixture. Enzymes specific to each of the four bacterial groups react with their specific chromogenic substrates to produce the colors that differentiate each bacterial type. In addition, the glucuronidase cleaves the fluorescent entity which results in a bright blue fluorescence from <i>E. coli</i> colonies lighted by long wave U.V. light.
Composition:	Contains nutrients and chromogenic substrates. pH 7.2 ± 0.2
Procedure:	Use standard methodology for product being tested. Incubation time and temperature; 24-48 hrs at 35°C. It is possible to take accurate readings at 24 hours incubation time. Colors will intensify and colony size will be larger at 36-48 hours. The fluorescence of <i>E. coli</i> can generally be detected by 12-15 hr. incubation. Fluorescence may diffuse throughout the plate after 24 hrs making quantification by fluorescence impossible after that initial time period.
Interpretation:	<i>E. coli</i> * colonies will have a dark blue/purple overall colony color and fluoresce a bright blue under long wave UV light. Non-fecal coliforms (<i>Enterobacter</i> , <i>Citrobacter</i> , <i>Klebsiella spp.</i> and others) do not fluoresce and are light blue/grey to purplish in color. <i>Salmonella spp.</i> do not fluoresce and appear as green/teal colonies. <i>Aeromonas spp.</i> do not fluoresce and are pink to very light pink. Gram-positive organisms are generally inhibited.
Storage:	Freeze
Shelf Life:	1 year

ECA Check and Easygel are registered trademarks of Micrology Laboratories LLC.

Micrology Laboratories, LLC., Goshen, IN 574.533.3351